

Genna

Made for end-users - reliable power, electric height adjustment & secure storage



 LiFePO4 battery

 Electric height adjustment

This cart supports laptops, micro PCs, and medical-grade all-in-one PCs. The configuration above shows a micro PC powering two monitors (up to 24") with independently adjustable monitor mounts to find the right amount of tilt for each user's ergonomic preference.

The dual screen setup lets users move quickly between apps, reducing time lost to switching, remembering details, or recovering after interruptions. These small delays add up, so minimising them improves workflow and frees up more time to care. It also enables multiple clinicians to review patient information at the same time.

For more info on our carts, please get in touch with our sales team, who will be happy to help. Our contact details are below.

Features

An electric height-adjustable cart with a generous storage capacity, a large work surface for medicines prep, and a unique, ergonomic design to enhance accessibility to patient information and IT hardware.

- ✓ Electric height adjustment to change height at the touch of a button, with independent height adjustment for your monitor.
- ✓ Secure, organised storage with up to five drawers reduces the need for restocking. Digi, keyed, and RFID lock options are available.
- ✓ Integrated bottle drawer allows for up to six 500ml bottles of liquid medication to be stored separately on two levels.
- ✓ Retractable keyboard conveniently tucks out of the way during medication prep and dispensing.
- ✓ Work surface is durable and resistant to hospital cleaning agents.

Benefits

- Small footprint makes the cart easy to position in busy bays.
- Large work surface and slide out tray provide ample medication preparation space.
- Conveniently placed grab handles at the side of the cart provide an unobscured view when moving the cart to a new area.
- The wider work surface design offers protection from damage for your screen while being moved to a new area.
- LED battery gauge at front of cart provides an at-a-glance view of battery power status from across the room.
- Easily accessible PC enclosure and cable management.
- Up to two optional side rails provide versatility for mounting a range of accessories, such as label printers, barcode scanners, and sharps bins.
- Designed for infection control and ease of cleaning.
- Comprehensive 5-year cart warranty.

Areas of use

Genna is ideal for providing access to electronic patient records as well as secure access to medication and supplies in areas such as:

- A&E / ED
- ICU
- Theatres
- Wards

Compatible with



Infinity
all-in-one PC



Micro
PC




Laptop




Dual
screens

Drawings



PC caddy tray location and dimensions



Main specifications

Work surface height (extended)	1283 mm
Work surface height (compressed)	883 mm
Work surface dimensions	646 mm (width) x 549 mm (depth)
Weight (cart base only)	76.6 kg
Weight (max weight possible)	98 kg
PC caddy dimensions	343 mm (width) x 234 mm (depth) x 101 mm (height)
PC shelf dimensions	196 mm (width) x 234 mm (depth) x 52 mm (height)

Approved cleaning solutions

- Quaternary ammonium chloride compounds (such as AirX109)
- Bleach-Household strength (1:10 minimum dilution strength)
- Ethyl or Isopropyl alcohol-70%
- Hydrogen Peroxide
- Phenolic disinfectants (such as Lysol, WexCide)
- Ovation
- Fantastik
- Formula 409
- 1.6% Aqueous Ammonia
- "Green Soap" USP (United States Pharmacopoeia)
- Clinell Wipes
- Chlor Clean



Direct purchase



Made in the UK



Project management



Device configuration



Phased delivery



Instant support

Common accessories



Printer



Scanner



Medical keyboard & mouse

