

# medDispense® C series

## Maximum flexibility in a single cabinet

### Overview

The medDispense® C series automated dispensing cabinets deliver maximum flexibility in a single cabinet—from a single item to bulk supply storage. Easily expand your storage capacity further by combining the C series with CS series, CT series, F series, or V series.

- ✓ Units allow you to choose from a bulk supply, drawers, or a combination of both
- ✓ Drawers come in linear or matrix, depending on security and capacity requirements
- ✓ Hold from 112 items up to 600 items in 12 to 60 drawers, depending on the configuration
- ✓ Powered by medLogic™ to deliver increased efficiency and flexibility for your team and a safer environment for your patients

### Features

- Provides customizable capacity with both single item control and bulk storage capabilities, all in one cabinet
- Integrates with CS, CT, F, and V series medDispense® automated dispensing cabinets
- Linear or matrix drawer options to choose from, based on security level and item size
- Microban® antimicrobial infused surfaces
- Ability to easily combine bulk items with bins, and multiple drawer options in one cabinet
- Our smallest unit fits on a countertop for areas with lower medication capacity needs
- Comprehensive 5- year warranty

### Recommended for

- Psychiatric
- Critical access hospitals
- Simulation labs
- Long-term acute care

### Other unit sizes in this series



### Compatible with



CS series



CT series



F series



V series

# medDispense® C series

## Maximum flexibility in a single cabinet

### Accessories



Remote lock



Remote lock with temperature monitoring



External return bin



Barcode scanner



Bulk storage bins




Thermal printer



Label printer



Fingerprint scanner



Return drawer



Overbridge

### Specifications

Max unit capacity	112, 600, 700, or 800
Width	686 mm
Depth	559 mm (112-capacity unit and all compact units) or 711 mm (600-, 700-, and 800-capacity full depth units)
Height	1,550 mm
Weight	66-73 kgs (112-capacity unit), 90-101 kgs (600-capacity unit), 159-222 kgs (700-capacity unit), or 221-344 kgs (800-capacity unit)
Max number of drawers	12, 30, or 60 (dependent on unit capacity)
Material construction	Steel frame with polymer drawers and doors
Power	120VAC @ 60Hz

### Software

- All medDispense® automated dispensing cabinets utilize the medLogic™ software platform
- The medLogic™ software platform delivers increased efficiency and flexibility for your team and a safer environment for your patients.
- **Efficiency:**
  - Single medication dispensing across multiple medDispense® automated dispensing cabinets.
  - Restock faster with barcode stocking.
  - FIFO stocking workflow ensures oldest stocked items are dispensed first, which helps reduced expired medication expenses.
- **Safety:**
  - Fully supports the UK NHS ECD register.
  - On-demand stock report generation by item/medication.
  - Help reduce medication errors with our configurable medication search.
  - Increase medication traceability with lot number stocking.
- **Flexibility:**
  - Fractional dose stocking & dispensing helps reduce medication waste.
  - High definition 16:9 and laptop display with the PharmacyCenter™ and NurSync™ modules provides increased technology flexibility for the pharmacy team.
- medLogic™ secures and tracks both your medications and supplies while seamlessly integrating with any hospital information system.
- medLogic™ modules improve workflow efficiency, safety and visibility so you can save time and focus on delivering exceptional patient care.
- **NurSync™ Remote Medication Queuing:**
  - Enables nurses to create medication queues for their patients to help reduce medication retrieval times, create a more efficient workflow and spend more time with patients.
- **medGuard™ Controlled Drug Inventory Management:**
  - Helps close the loop of controlled drug medication delivery and monitors for potential diversion with controlled restocking access and software safeguards.
  - Discrepancy tracking ensures the retrieval quantity recorded in medGuard™ matches the quantity restocked in the medDispense® automated dispensing cabinet.
- **PharmacyCenter™ Reporting and System Administration:**
  - Enables pharmacy staff to manage their medications across multiple medDispense® automated dispensing cabinets from one easy-to-use dashboard and reporting system.
  - Manage and maintain your facility's automated dispensing cabinets from any PC.
  - Visibility and management of your inventory across all automated dispensing cabinets.
  - Report on key metrics that matter to your facility, such as: inventory, top usage, transactions, controlled medications, and more.

