

SurgiView

Mobilise large-format surgical displays in theatres



Designed for Surgeons, showcased at the Future Surgery Show



IN PARTNERSHIP WITH



Royal College
of Surgeons
of England
ADVANCING SURGICAL CARE

Features

The SurgiView is designed to mobilise large screen surgical displays in theatres. This flexible approach is particularly beneficial in hybrid Operating Rooms and for use in Robotic Surgery Theatres.

We frequently see challenges with mounting on/in-wall solutions in just the right position, especially in PFI buildings. The SurgiView addresses mounting and power restraints and reduces downtime for installation, providing your clinical teams with the versatility they need.

- ✓ Compatible with screens up to 55"
- ✓ Height-adjustable keyboard tray
- ✓ Secure storage for PCs
- ✓ Capable of housing a medical grade isolation transformer, or our EN 60601 compliant power supply for mobile use-cases.
- ✓ Made from robust and durable materials that are resistant to hospital cleaning agents


Benefits

- Large, lockable castors provide high maneuverability.
- Can be supplied with a compliant display, ensuring compliance to RCR PACs clinical review standards, or other relevant guidelines.
- Options for medical keyboards, mice and other peripherals.
- With Parity Medical's configuration service, the cart can be fitted with trust equipment and neatly cable managed, with screen installation on-site. This service both reduces strain on hospital resources, whilst speeding up the roll-out process.
- Designed for infection control and ease of cleaning.
- Comprehensive 5-year cart warranty.

Areas of use

- Theatres
- MDT
- Teaching environments
- Patient status at-a-glance/bed management


Recommended displays



Parity Medical
MedView



Barco surgical
monitors




Eizo Curator
range

SurgiView

Mobilise large-format surgical displays in theatres

Drawings



Main specifications

Rear work surface height 973 mm

Cart base dimensions 620 mm (width) x 774 mm (depth)

Approved cleaning solutions

- Quaternary ammonium chloride compounds (such as AirX109)
- Bleach-Household strength (1:10 minimum dilution strength)
- Ethyl or Isopropyl alcohol-70%
- Hydrogen Peroxide
- Phenolic disinfectants (such as Lysol, WexCide)
- Ovation
- Fantastic
- Formula 409
- 1.6% Aqueous Ammonia
- "Green Soap" USP (United States Pharmacopoeia)
- Clinell Wipes
- Chlor Clean



Direct purchase



Made in the UK



Project management



Device configuration



Phased delivery



Instant support

Common accessories



Medical keyboard & mouse



File storage

Version 1.0