

# MedStore Plus

## The big one! Our largest-capacity workstation

### Features

Carry more medication with the MedStore Plus! With increased storage capacity, wider drawers, and a large work surface, this cart provides an ideal central workstation for ePMA at the point of care. A large number of drawers make this solution great for organising medication and reducing the need to restock, saving you time while working.

- ✓ High-storage capacity workstation for medication dispensing at the point of care.
- ✓ Can be configured to allow for the addition of a range of accessories and drawer configurations.
- ✓ Work surface is durable and resistant to hospital cleaning agents.

### Benefits

- Large work surface for preparing drugs and paper documents.
- Slide-out mouse tray for left or right-handed use helps to save space on the work surface.
- Conveniently placed grab handles allow precise and easy movement.
- Castors provide smooth and stable movement even across uneven surfaces, such as thresholds.
- Front two wheels are lockable and ensure the cart remains stationary when in use.
- Designed for infection control and ease of cleaning.
- Drawers and trays can be configured to different depths.
- Secure lockable storage allows for access to medication. Digi, keyed, and RFID lock options are available.
- Magnetic door attaches to the side of the cart when open for increased safety.
- Comprehensive 5-year cart warranty.

### Areas of use

The MedStore Plus cart is ideal for providing access to electronic patient records as well as secure access to large amounts of medication in areas such as:

- A&E / ED
- ICU
- Theatres
- Wards



This cart is designed to be used with different computers and devices. Our medical-grade Infinity all-in-one PC is shown on the cart above and features:

- A large, intuitive touchscreen.
- Dual hot-swappable batteries for 24/7 constant power.
- Powerful CPU, ultra-fast SSD and up to 64GB system memory.
- Unique battery management software with remote monitoring.
- Front panel is IP65 rated for protection from water and dust.
- 3-year warranty.

For more info on our Infinity all-in-one PCs or to enquire about a cart bundle package, please get in touch with our sales team, who will be happy to help. Our contact details are below.

### Compatible with



Infinity all-in-one PC



Micro PC\*




Laptop

\*Requires SmartPower

# MedStore Plus

The big one! Our largest-capacity workstation

## Drawings



## Main specifications

Work surface height	1048 mm
Weight (base cart only)	62 kg
Weight (base cart with SmartPower)	76 kg
Work surface dimensions	669 mm (width) x 635 mm (depth)

## Approved cleaning solutions

- Quaternary ammonium chloride compounds (such as AirX109)
- Bleach-Household strength (1:10 minimum dilution strength)
- Ethyl or Isopropyl alcohol-70%
- Hydrogen Peroxide
- Phenolic disinfectants (such as Lysol, WexCide)
- Ovation
- Fantastic
- Formula 409
- 1.6% Aqueous Ammonia
- "Green Soap" USP (United States Pharmacopoeia)
- Clinell Wipes
- Chlor Clean



Direct purchase



Made in the UK



Project management



Device configuration



Phased delivery



Instant support

## Common accessories



Printer



Scanner



Medical keyboard & mouse

Version 2.1

