

# VersaLite EPR

## Ergonomic, durable, and built to last



This cart is designed to be used with different computers and devices. Our medical-grade Infinity all-in-one PC is shown on the cart above and features:

- A large, intuitive touchscreen.
- Dual hot-swappable batteries for 24/7 constant power.
- Powerful CPU, ultra-fast SSD and up to 64GB system memory.
- Unique battery management software with remote monitoring.
- Front panel is IP65 rated for protection from water and dust.
- 3-year warranty.

For more info on our Infinity all-in-one PCs or to enquire about a cart bundle package, please get in touch with our sales team, who will be happy to help. Our contact details are below.

## Features

The ergonomic nurse-led design of this cart provides the ultimate combination of functionality, mobility, and user comfort. It allows you to release more time to care, with instant access to patient records and a reduced risk of repetitive strain at the point of care.

- ✓ Simple hand-operated height adjustment.
- ✓ Retractable comfort keyboard tray with adjustable height, distance, and angle.
- ✓ Modular design allows for the addition of a range of accessories and drawer configurations.
- ✓ Work surface is durable and resistant to hospital cleaning agents.

## Benefits

- Small footprint makes the cart easy to position in busy bays.
- Conveniently placed grab handles allow precise and easy movement.
- Large rear castors provide smooth and stable movement even across uneven surfaces, such as thresholds.
- Designed for infection control and ease of cleaning.
- An ergonomically designed base allows the user to find their perfect standing or seated position.
- Comfort keyboard tray is retractable and features a simple lift and release mechanism allowing the user to set the height, distance, and angle independently of the workstation, so the user can maintain a healthy posture.
- Comprehensive 5-year cart warranty.

## Areas of use

This EPR cart is ideal for providing access to electronic patient records or medication notes in areas such as:

- A&E / ED
- ICU
- Theatres
- Wards

## Compatible with



Infinity  
all-in-one PC




Laptop

# VersaLite EPR

Ergonomic, durable, and built to last

## Drawings



## Main specifications

Work surface height (extended)	1218 mm
Work surface height (compressed)	864 mm
Weight (base cart only)	34.9 kg
Work surface dimensions	513 mm (width) x 636 mm (depth)
Cart base dimensions	515 mm (width) x 627 mm (depth)

## Approved cleaning solutions

- Quaternary ammonium chloride compounds (such as AirX109)
- Bleach-Household strength (1:10 minimum dilution strength)
- Ethyl or Isopropyl alcohol-70%
- Hydrogen Peroxide
- Phenolic disinfectants (such as Lysol, WexCide)
- Ovation
- Fantastic
- Formula 409
- 1.6% Aqueous Ammonia
- "Green Soap" USP (United States Pharmacopoeia)
- Clinell Wipes
- Chlor Clean



Direct purchase



Made in the UK



Project management



Device configuration



Phased delivery



Instant support

## Common accessories



Printer



Scanner



Medical keyboard & mouse

Version 2.1

